## Siebreinigungsmodul SRH 4/200

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## Sieve cleaning module SRH 4/200

### Intended use

Device for the cleaning of laboratory and analysing sieves in an ultrasonic cleaning unit\*.

\*Suitable Elma ultrasonic cleaning units:  
Elmasonic S300/H; S450H; S900H.

### User

Instructed staff only

### Safety instructions

	<ul style="list-style-type: none"><li>• Risk of crushing injury due to the top part falling down!</li><li>• Transportation of the sieve cleaning module and placement of the module into the ultrasonic tank or removal out of the ultrasonic tank with closed top part only!</li></ul>
	<ul style="list-style-type: none"><li>• Risk of injury during operation of the sieve cleaning module due to the rotating shaft inside! Hair or loose clothing may be drawn in.</li><li>• Operation of the sieve cleaning module with closed top part only!</li><li>• Switch off the sieve cleaning module before you open the top part of the sieve cleaning module!</li></ul>

### Technical Details

Mains voltage (V~)	100-120 oder 220-240
Voltage supply drive (V~)	24
Power consumption (W)	10
Overall dimensions W/D/H (mm)	580/260/300
Weight (kg)	ca. 6,8

### How to proceed

1. Fill the ultrasonic tank up to the marking on the tank inside (filling level of approx 15 cm in an Elmasonic S 300) with water and suitable cleaning medium. Operate the ultrasound of approx 10 min to degas the cleaning liquid (please see also Operating Instructions of the ultrasonic unit). Heat up the cleaning liquid if required.

2. Place the sieves into the sieve cleaning module and close the top part of the module. Insert the sieve cleaning module into the tank of the ultrasonic unit. Cleaning of up to 4 sieves [Ø max. 200/203 mm (8")] in one cleaning cycle. Sieves of different brands can be cleaned simultaneously.

3. Connect the mains cable to the mains supply socket.  
Switch on the rotation of the sieves at the rocker switch.  
Switch on the ultrasound.  
Check the cleaning result by visual inspection (5 – 10 minutes depending on grain properties), repeat the cleaning cycle if required.  
Switch off the ultrasonic unit and the sieve cleaning module after termination of the cleaning cycle.

4. Take out individual sieves or the complete sieve cleaning module.  
The plastic cover of the ultrasonic unit can be used as a drip off collecting tray.

The sieve cleaning module is in compliance with the CE marking criteria with regard to the Directives 2006 / 95 / EWG (EC / EEC) and 2004 / 108 / EWG (EC / EEC).

### Putting out of action and waste disposal

The unit components can be taken to metal and electronics recycling stations or returned to the manufacturer.

### Manufacturer's contact address

**Elma Schmidbauer GmbH**  
Gottlieb-Daimler-Str. 17  
78224 Singen (Deutschland)  
[www.elma-ultrasonic.com](http://www.elma-ultrasonic.com)

### Technical Support

Tel: +49 (0) 77 31 / 882-280  
E-Mail: [support@elma-ultrasonic.com](mailto:support@elma-ultrasonic.com)

Copyright © 2020 Elma Schmidbauer GmbH.  
All rights reserved.  
Subject to technical and visual modifications.

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